



Mini Hippocampal Tool - 14.5cm



Both ends are malleable to allow you to adjust their shape to the correct angle for optimal positioning.

General Info

Product Number	10099-12
Category	Spoons & Spatulae
EAN	4251909500788



Find 10099-12 on
finescience.com

Technical Info

Shipping Weight	25 g	Net Weight	22 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Length	14.5 cm		
Tipdimensions	20 x 5 mm & 13 x 4 mm		
Tipthickness	0.5 mm		