



Diamond Knife - Angled 30°/12cm/1mm Cutting Edge



This knife has a gem-quality, neutral diamond blade which is honed to the finest cutting edge possible. The blade is made to an ophthalmic surgical standard and is inspected under 400x magnification. The blade retracts completely down into the handle for safe care and storage. Cautionary note: This blade is extremely fragile and vulnerable to lateral pressure.

General Info

Product Number	10100-30
Category	Scalpels & Knives > Micro Knives
EAN	4251909500801



Find 10100-30 on finescience.com

Technical Info

Shipping Weight	55 g	Net Weight	10 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Anodized aluminum handle with diamond blade		
Antimagnetic	Yes		
Autoclavesafe	Yes		
Bladeshape	Lancet		
Cuttingangle	30 °		
Cuttingedge	1 mm		
Length	12 cm		
Tipthickness	0.15 mm		