



Tungsten Needles - 0.25mm Rod Diameter



For very fine micro dissection, providing lateral stability as well as minimal spring action during manipulation.

General Info

Product Number	10130-10
Category	Pins, Probes & Hooks > Pins Most Popular > Pins, Probes & Hooks
EAN	4251909500870



Find 10130-10 on
finescience.com

Technical Info

Shipping Weight	30 g	Net Weight	5 g
Package Size	10 Pc.	Unit of Measure	PA
Minimum Order Quantity	1		
Alloymaterial	Tungsten Carbide		
Length	1.2 cm		
Roddiameter	0.25 mm		
Tipdiameter	1 µm		