



Probe - Angled 80° Short/15.5cm/0.15mm Tip Diameter



General Info

Product Number	10140-02
Category	Pins, Probes & Hooks > Probes Most Popular > Pins, Probes & Hooks
EAN	4061706158873



Find 10140-02 on
finescience.com

Technical Info

Shipping Weight	20 g	Net Weight	15 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Length	15.5 cm		
Tipdiameter	0.15 mm		
Tipshape	Angled 80°		