



Probe - Angled 80° Long/15.5cm/0.15mm Tip Diameter



Allgemein

Produktnummer	10140-03
Kategorie	Pins, Probes & Hooks > Probes Most Popular > Pins, Probes & Hooks
EAN	4251909500894



Finde 10140-03 auf
finescience.com

Technische Informationen

Versandgewicht	20 g	Nettogewicht	15 g
Paketgröße	1 Pc.	Maßeinheit	EA
Mindestbestellmenge	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Length	15.5 cm		
Tipdiameter	0.15 mm		
Tipshape	Angled 80° Long		