



Probe - Angled 80° Long/15.5cm/0.15mm Tip Diameter



General Info

Product Number	10140-03
Category	Pins, Probes & Hooks > Probes Most Popular > Pins, Probes & Hooks
EAN	4251909500894



Find 10140-03 on
finescience.com

Technical Info

Shipping Weight	20 g	Net Weight	15 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Length	15.5 cm		
Tipdiameter	0.15 mm		
Tipshape	Angled 80° Long		