



Gross Anatomy Knife - Straight/30cm/15cm Cutting Edge



Made of molybdenum stainless steel to keep the cutting edge sharp. Molybdeum is an alloying agent used to form extremely hard steel with greater corrosion resistance. The yellow handles are fiberglass reinforced nylon with slip resistant surface.

General Info

Product Number	10157-15
Category	Scalpels & Knives > Knives
EAN	4033182001660



Find 10157-15 on finescience.com

Technical Info

Shipping Weight	130 g	Net Weight	122 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Stainless Steel (Inox)		
Autoclavesafe	No		
Bladeshape	Straight		
Cuttingedge	15 cm		
Length	28 cm		
Thickness	2.15 mm		