



Titanium Y-Hook - Angled 45°/12.5cm



Lightweight, color coded for easy identification, have an ergonomic design and a non-reflective surface.

General Info

Product Number	10166-11
Category	Pins, Probes & Hooks > Hooks
EAN	4251909500948



Find 10166-11 on
finescience.com

Technical Info

Shipping Weight	30 g	Net Weight	9 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Titanium		
Antimagnetic	Yes		
Length	12.5 cm		
Tipthickness	0.5 mm		
Tipwidth	1 mm		
Tipshape	Angled 45°		