



Grooved Director with Tongue, 13 cm



The Grooved Director is a specialized multipurpose device which is used by surgeons in different areas to probe cavities and wounds, enlarge fascia and isolate the tongue for example.

General Info

Product Number	10290-13
Category	Pins, Probes & Hooks > Probes
EAN	4251909511531



Find 10290-13 on
finescience.com

Technical Info

Shipping Weight	7 g	Net Weight	6 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Length	13 cm		
Sterile	No		