



Micro Knife - Plastic Handle/Crescent Blade/14cm



MOST
POPULAR



NOT
AUTOCLAVABLE

Thin blades make these ideal for the finest incisions. The 12.5cm plastic handle has a slip-proof grip. Crescent knives are designed to create a smooth scleral shelf to promote a self-sealing wound. The blade’s consistent sharpness virtually eliminates tissue drag and distortion.

General Info

Product Number	10317-14
Category	Scalpels & Knives > Micro Knives Most Popular > Scalpels & Knives
EAN	4251909500986



Find 10317-14 on
finescience.com

Technical Info

Shipping Weight	10 g	Net Weight	2 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Carbon steel tip / Plastic handle		
Autoclavesafe	No		
Bladeshape	Crescent blade		
Cuttingedge	11 mm		
Length	14 cm		
Sterile	Yes		
Tipthickness	0.2 mm		