



Micro Knife - Plastic Handle/Needle Blade/14.5cm



MOST
POPULAR



NOT
AUTOCLAVABLE

Thin blades make these ideal for the finest incisions. The 12.5cm plastic handle has a slip-proof grip. Specialty knives are available in a variety of shapes and sizes for various techniques and procedures.

General Info

Product Number	10318-14
Category	Scalpels & Knives > Micro Knives Most Popular > Scalpels & Knives
EAN	4251909500993



Find 10318-14 on
finescience.com

Technical Info

Shipping Weight	15 g	Net Weight	2 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Carbon steel tip / Plastic handle		
Autoclavesafe	No		
Bladeshape	Needle blade		
Cuttingedge	4 mm		
Length	14.5 cm		
Sterile	Yes		
Tipthickness	0.2 mm		