



Mirror Finish Cutting Spoon - 13cm/4.5mm Tip Diameter



Non stick cutting spoon.

General Info

Product Number	10360-13
Category	Spoons & Spatulae
EAN	4251909501013



Find 10360-13 on
finescience.com

Technical Info

Shipping Weight	20 g	Net Weight	16 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Length	13 cm		
Spoonddepth	0.15 mm		
Tipdiameter	4.5 mm		