



Graefe Iris Knife - Straight/13cm



General Info

Product Number	10370-13
Category	Scalpels & Knives > Micro Knives
EAN	4251909501020



Find 10370-13 on
finescience.com

Technical Info

Shipping Weight	50 g	Net Weight	15 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Autoclavesafe	Yes		
Bladeshape	Straight		
Cuttingedge	28 mm		
Length	13 cm		
Thickness	0.6 mm		