



Moria MC 17 BIS Perforated Spoon - 14.5cm/15mm Tip Diameter



Most commonly used to take embryos out of solution and handling very small organs.

General Info

Product Number	10370-18
Category	Spoons & Spatulae Most Popular > Spoons & Spatulae
EAN	4251909501044



Find 10370-18 on
finescience.com

Technical Info

Shipping Weight	20 g	Net Weight	17 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Autoclavesafe	Yes		
Length	14.5 cm		
Holediameter	0.9 mm		
Spoonddepth	3 mm		
Tipdiameter	15 mm		