



Moria MC17C Mini Perforated Spoon - 10cm/8mm Tip Diameter



Most commonly used to take embryos out of solution and handling very small organs.

General Info

Product Number	10370-19
Category	Spoons & Spatulae
EAN	4251909501051



Find 10370-19 on
finescience.com

Technical Info

Shipping Weight	15 g	Net Weight	7 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Length	10 cm		
Holediameter	0.8 mm		
Spoonddepth	1 mm		
Tipdiameter	8 mm		