



Anatomical Standard Forceps - Serrated/Straight/12cm



Allgemein

Produktnummer	11000-12
Kategorie	Forceps & Tweezers > Serrated
EAN	4061706023065



Finde 11000-12 auf finescience.com

Technische Informationen

Versandgewicht	20 g	Nettogewicht	18 g
Paketgröße	1 Pc.	Maßeinheit	EA
Mindestbestellmenge	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Length	12 cm		
Tipdimensions	2.5 x 1.35 mm		
Tipwidth	2.5 mm		
Tips	Serrated		
Tipshape	Straight		