



Anatomical Standard Forceps - Serrated/Straight/13cm



Allgemein

Produktnummer	11000-13
Kategorie	Forceps & Tweezers > Serrated
EAN	4251909501068



Finde 11000-13 auf finescience.com

Technische Informationen

Versandgewicht	25 g	Nettogewicht	20 g
Paketgröße	1 Pc.	Maßeinheit	EA
Mindestbestellmenge	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Autoclavesafe	Yes		
Length	13 cm		
Tipdimensions	2.8 x 1.4 mm		
Tipwidth	2.8 mm		
Tips	Serrated		
Tipshape	Straight		