



Anatomical Standard Forceps - Serrated/Straight/14.5cm



General Info

Product Number	11000-14
Category	Forceps & Tweezers > Serrated
EAN	4251909501075



Find 11000-14 on
finescience.com

Technical Info

Shipping Weight	30 g	Net Weight	28 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Length	14.5 cm		
Tipdimensions	3.25 x 1.6 mm		
Tipwidth	3.25 mm		
Tips	Serrated		
Tipshape	Straight		