



Anatomical Standard Forceps - Serrated/Straight/16cm



Allgemein

Produktnummer	11000-16
Kategorie	Forceps & Tweezers > Serrated
EAN	4061706023126



Finde 11000-16 auf finescience.com

Technische Informationen

Versandgewicht	35 g	Nettogewicht	32 g
Paketgröße	1 Pc.	Maßeinheit	EA
Mindestbestellmenge	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Autoclavesafe	Yes		
Length	16 cm		
Tipdimensions	3.3 x 1.7 mm		
Tipwidth	3.3 mm		
Tips	Serrated		
Tipshape	Straight		