



Anatomical Standard Forceps - Serrated/Straight/16cm



General Info

Product Number	11000-16
Category	Forceps & Tweezers > Serrated
EAN	4061706023126



Find 11000-16 on
finescience.com

Technical Info

Shipping Weight	35 g	Net Weight	32 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Length	16 cm		
Tipdimensions	3.3 x 1.7 mm		
Tipwidth	3.3 mm		
Tips	Serrated		
Tipshape	Straight		