



Anatomical Standard Forceps - Serrated/Straight/20cm



Allgemein

Produktnummer	11000-20
Kategorie	Forceps & Tweezers > Serrated
EAN	4061706023140



Finde 11000-20 auf finescience.com

Technische Informationen

Versandgewicht	55 g	Nettogewicht	50 g
Paketgröße	1 Pc.	Maßeinheit	EA
Mindestbestellmenge	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Length	20 cm		
Tipdimensions	3.75 x 1.9 mm		
Tipwidth	3.75 mm		
Tips	Serrated		
Tipshape	Straight		