



Anatomical Standard Forceps - Serrated/Straight/30cm



Allgemein

Produktnummer	11000-30
Kategorie	Forceps & Tweezers > Serrated
EAN	4061706023188



Finde 11000-30 auf finescience.com

Technische Informationen

Versandgewicht	105 g	Nettogewicht	98 g
Paketgröße	1 Pc.	Maßeinheit	EA
Mindestbestellmenge	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Length	30 cm		
Tipdimensions	4 x 2.25 mm		
Tipwidth	4 mm		
Tips	Serrated		
Tipshape	Straight		