



Anatomical Standard Forceps - Serrated/Straight/30cm



General Info

Product Number	11000-30
Category	Forceps & Tweezers > Serrated
EAN	4061706023188



Find 11000-30 on
finescience.com

Technical Info

Shipping Weight	105 g	Net Weight	98 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Length	30 cm		
Tipdimensions	4 x 2.25 mm		
Tipwidth	4 mm		
Tips	Serrated		
Tipshape	Straight		