



Anatomical Standard Forceps - Serrated/Curved/13cm



Allgemein

Produktnummer	11001-13
Kategorie	Forceps & Tweezers > Serrated
EAN	4251909501099



Finde 11001-13 auf finescience.com

Technische Informationen

Versandgewicht	20 g	Nettogewicht	18 g
Paketgröße	1 Pc.	Maßeinheit	EA
Mindestbestellmenge	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Length	13 cm		
Tipdimensions	2.8 x 1.4 mm		
Tipwidth	2.8 mm		
Tips	Serrated		
Tipshape	Curved		