



Anatomical Standard Forceps - Serrated/Curved/13cm



General Info

Product Number	11001-13
Category	Forceps & Tweezers > Serrated
EAN	4251909501099



Find 11001-13 on
finescience.com

Technical Info

Shipping Weight	20 g	Net Weight	18 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Length	13 cm		
Tipdimensions	2.8 x 1.4 mm		
Tipwidth	2.8 mm		
Tips	Serrated		
Tipshape	Curved		