



Anatomical Standard Forceps - Serrated/Curved/20cm



General Info

|                |                               |
|----------------|-------------------------------|
| Product Number | 11001-20                      |
| Category       | Forceps & Tweezers > Serrated |
| EAN            | 4061706023287                 |



Find 11001-20 on [finescience.com](https://finescience.com)

Technical Info

|                        |                                 |                 |      |
|------------------------|---------------------------------|-----------------|------|
| Shipping Weight        | 50 g                            | Net Weight      | 46 g |
| Package Size           | 1 Pc.                           | Unit of Measure | EA   |
| Minimum Order Quantity | 1                               |                 |      |
| Alloymaterial          | Surgical Stainless Steel (Inox) |                 |      |
| Length                 | 20 cm                           |                 |      |
| Tipdimensions          | 3.75 x 1.9 mm                   |                 |      |
| Tipwidth               | 3.75 mm                         |                 |      |
| Tips                   | Serrated                        |                 |      |
| Tipshape               | Curved                          |                 |      |