



Anatomical Narrow Forceps - Serrated/Straight/14.5cm



Lighter spring action than standard pattern for grasping more delicate tissue.

General Info

Product Number	11002-14
Category	Forceps & Tweezers > Serrated Most Popular > Forceps & Tweezers
EAN	4251909501143



Find 11002-14 on
finescience.com

Technical Info

Shipping Weight	30 g	Net Weight	25 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Length	14.5 cm		
Tipdimensions	2.25 x 1.45 mm		
Tipwidth	2.25 mm		
Tips	Serrated		
Tipshape	Straight		