



Anatomical Narrow Forceps - Serrated/Straight/18cm



Lighter spring action than standard pattern for grasping more delicate tissue.

General Info

Product Number	11002-18
Category	Forceps & Tweezers > Serrated Most Popular > Forceps & Tweezers
EAN	4251909501167



Find 11002-18 on
finescience.com

Technical Info

Shipping Weight	40 g	Net Weight	38 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Length	18 cm		
Tipdimensions	2.4 x 1.6 mm		
Tipwidth	2.4 mm		
Tips	Serrated		
Tipshape	Straight		