



Anatomical Narrow Forceps - Serrated/Straight/25cm



Lighter spring action than standard pattern for grasping more delicate tissue.

General Info

Product Number	11002-25
Category	Forceps & Tweezers > Serrated Most Popular > Forceps & Tweezers
EAN	4251909501181



Find 11002-25 on
finescience.com

Technical Info

Shipping Weight	70 g	Net Weight	67 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Length	25 cm		
Tipdimensions	3.2 x 2 mm		
Tipwidth	3.2 mm		
Tips	Serrated		
Tipshape	Straight		