



Anatomical Narrow Forceps - Serrated/Curved/12cm



Lighter spring action than standard pattern for grasping more delicate tissue.

General Info

|                |                                                                    |
|----------------|--------------------------------------------------------------------|
| Product Number | 11003-12                                                           |
| Category       | Forceps & Tweezers > Serrated<br>Most Popular > Forceps & Tweezers |
| EAN            | 4251909501198                                                      |



Find 11003-12 on  
[finescience.com](https://finescience.com)

Technical Info

|                        |                                 |                 |      |
|------------------------|---------------------------------|-----------------|------|
| Shipping Weight        | 15 g                            | Net Weight      | 14 g |
| Package Size           | 1 Pc.                           | Unit of Measure | EA   |
| Minimum Order Quantity | 1                               |                 |      |
| Alloymaterial          | Surgical Stainless Steel (Inox) |                 |      |
| Length                 | 12 cm                           |                 |      |
| Tipdimensions          | 2 x 1.25 mm                     |                 |      |
| Tipwidth               | 2 mm                            |                 |      |
| Tips                   | Serrated                        |                 |      |
| Tipshape               | Curved                          |                 |      |