



Anatomical Narrow Forceps - Serrated/Curved/13cm



Lighter spring action than standard pattern for grasping more delicate tissue.

General Info

Product Number	11003-13
Category	Forceps & Tweezers > Serrated Most Popular > Forceps & Tweezers
EAN	4251909501204



Find 11003-13 on
finescience.com

Technical Info

Shipping Weight	20 g	Net Weight	18 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Length	13 cm		
Tipdimensions	2.15 x 1.3 mm		
Tipwidth	2.15 mm		
Tips	Serrated		
Tipshape	Curved		