



Anatomical Narrow Forceps - Serrated/Curved/16cm



Lighter spring action than standard pattern for grasping more delicate tissue.

General Info

Product Number	11003-16
Category	Forceps & Tweezers > Serrated Most Popular > Forceps & Tweezers
EAN	4251909501228



Find 11003-16 on
finescience.com

Technical Info

Shipping Weight	30 g	Net Weight	29 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Autoclavesafe	Yes		
Length	16 cm		
Tipdimensions	2.35 x 1.5 mm		
Tipwidth	2.35 mm		
Tips	Serrated		
Tipshape	Curved		