



Anatomical Narrow Forceps - Serrated/Curved/18cm



Lighter spring action than standard pattern for grasping more delicate tissue.

General Info

Product Number	11003-18
Category	Forceps & Tweezers > Serrated Most Popular > Forceps & Tweezers
EAN	4251909501235



Find 11003-18 on
finescience.com

Technical Info

Shipping Weight	40 g	Net Weight	36 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Length	18 cm		
Tipdimensions	2.4 x 1.6 mm		
Tipwidth	2.4 mm		
Tips	Serrated		
Tipshape	Curved		