



Anatomical Narrow Forceps - Serrated/Curved/20cm



Lighter spring action than standard pattern for grasping more delicate tissue.

General Info

Product Number	11003-20
Category	Forceps & Tweezers > Serrated Most Popular > Forceps & Tweezers
EAN	4251909501242



Find 11003-20 on
finescience.com

Technical Info

Shipping Weight	45 g	Net Weight	45 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Length	20 cm		
Tipdimensions	2.75 x 1.8 mm		
Tipwidth	2.75 mm		
Tips	Serrated		
Tipshape	Curved		