



Adson Forceps - Serrated/Straight/12cm



General Info

Product Number	11006-12
Category	Forceps & Tweezers > Serrated Most Popular > Forceps & Tweezers
EAN	4251909501259



Find 11006-12 on
finescience.com

Technical Info

Shipping Weight	25 g	Net Weight	18 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Length	12 cm		
Tipdimensions	1.5 x 1 mm		
Tipwidth	1.5 mm		
Tips	Serrated		
Tipshape	Straight		