



Vial Forceps- 15cm



Tip tapers for 6 cm, allowing easy access to specimens in vials.

General Info

Product Number	11006-15
Category	Forceps & Tweezers > Smooth
EAN	4251909501266



Find 11006-15 on finescience.com

Technical Info

Shipping Weight	50 g	Net Weight	17 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Length	15 cm		
Tipdiameter	0.5 mm		
Tipdimensions	0.5 x 0.4 mm		
Tips	Smooth		
Tipshape	Straight		