



Semken Anatomical Forceps - Serrated/Straight/13cm



Allgemein

Produktnummer	11008-13
Kategorie	Forceps & Tweezers > Serrated
EAN	4251909501273



Finde 11008-13 auf finescience.com

Technische Informationen

Versandgewicht	15 g	Nettogewicht	11 g
Paketgröße	1 Pc.	Maßeinheit	EA
Mindestbestellmenge	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Length	13 cm		
Tipdimensions	1.3 x 1 mm		
Tipwidth	1.3 mm		
Tips	Serrated		
Tipshape	Straight		