



Semken Anatomical Forceps - Serrated/Straight/13cm



General Info

Product Number	11008-13
Category	Forceps & Tweezers > Serrated
EAN	4251909501273



Find 11008-13 on
finescience.com

Technical Info

Shipping Weight	15 g	Net Weight	11 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Length	13 cm		
Tipdimensions	1.3 x 1 mm		
Tipwidth	1.3 mm		
Tips	Serrated		
Tipshape	Straight		