



Semken Anatomical Forceps - Serrated/Straight/15cm



Allgemein

Produktnummer	11008-15
Kategorie	Forceps & Tweezers > Serrated
EAN	4251909501280



Finde 11008-15 auf finescience.com

Technische Informationen

Versandgewicht	25 g	Nettogewicht	18 g
Paketgröße	1 Pc.	Maßeinheit	EA
Mindestbestellmenge	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Length	15 cm		
Tipdimensions	1.5 x 1.1 mm		
Tipwidth	1.5 mm		
Tips	Serrated		
Tipshape	Straight		