



Semken Anatomical Forceps - Serrated/Straight/15cm



General Info

Product Number	11008-15
Category	Forceps & Tweezers > Serrated
EAN	4251909501280



Find 11008-15 on
finescience.com

Technical Info

Shipping Weight	25 g	Net Weight	18 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Length	15 cm		
Tipdimensions	1.5 x 1.1 mm		
Tipwidth	1.5 mm		
Tips	Serrated		
Tipshape	Straight		