



Semken Anatomical Forceps - Serrated/Curved/13cm



Allgemein

|               |                               |
|---------------|-------------------------------|
| Produktnummer | 11009-13                      |
| Kategorie     | Forceps & Tweezers > Serrated |
| EAN           | 4251909501297                 |



Finde 11009-13 auf [finescience.com](https://finescience.com)

Technische Informationen

|                     |                                 |              |      |
|---------------------|---------------------------------|--------------|------|
| Versandgewicht      | 15 g                            | Nettogewicht | 11 g |
| Paketgröße          | 1 Pc.                           | Maßeinheit   | EA   |
| Mindestbestellmenge | 1                               |              |      |
| Alloymaterial       | Surgical Stainless Steel (Inox) |              |      |
| Length              | 13 cm                           |              |      |
| Tipdimensions       | 1.3 x 1 mm                      |              |      |
| Tipwidth            | 1.3 mm                          |              |      |
| Tips                | Serrated                        |              |      |
| Tipshape            | Curved                          |              |      |