



Taylor Forceps - Serrated/Curved/17cm



General Info

Product Number	11011-17
Category	Forceps & Tweezers > Serrated
EAN	4251909501327



Find 11011-17 on
finescience.com

Technical Info

Shipping Weight	25 g	Net Weight	24 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Length	17 cm		
Tipdimensions	1.5 x 1 mm		
Tips	Serrated		
Tipshape	Curved		