



Potts-Smith Forceps - Serrated/Straight/18cm



General Info

|                |                               |
|----------------|-------------------------------|
| Product Number | 11012-18                      |
| Category       | Forceps & Tweezers > Serrated |
| EAN            | 4251909501334                 |



Find 11012-18 on  
[finescience.com](https://www.finescience.com)

Technical Info

|                        |                                 |                 |      |
|------------------------|---------------------------------|-----------------|------|
| Shipping Weight        | 40 g                            | Net Weight      | 36 g |
| Package Size           | 1 Pc.                           | Unit of Measure | EA   |
| Minimum Order Quantity | 1                               |                 |      |
| Alloymaterial          | Surgical Stainless Steel (Inox) |                 |      |
| Length                 | 18 cm                           |                 |      |
| Tipdimensions          | 1.25 x 1 mm                     |                 |      |
| Tips                   | Serrated                        |                 |      |
| Tipshape               | Straight                        |                 |      |