



Potts-Smith Forceps - Serrated/Straight/21cm



General Info

Product Number	11012-21
Category	Forceps & Tweezers > Serrated
EAN	4251909501341



Find 11012-21 on
[finescience.com](https://www.finescience.com)

Technical Info

Shipping Weight	45 g	Net Weight	40 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Length	21 cm		
Tipdimensions	1.8 x 1.5 mm		
Tips	Serrated		
Tipshape	Straight		