



Potts-Smith Forceps - Serrated/Curved/21cm



General Info

Product Number	11013-21
Category	Forceps & Tweezers > Serrated
EAN	4251909501372



Find 11013-21 on
finescience.com

Technical Info

Shipping Weight	45 g	Net Weight	40 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Length	21 cm		
Tipdimensions	1.8 x 1.5 mm		
Tips	Serrated		
Tipshape	Curved		