



Gerald Forceps - Serrated/Straight/18cm



General Info

Product Number	11014-18
Category	Forceps & Tweezers > Serrated
EAN	4251909501389



Find 11014-18 on
finescience.com

Technical Info

Shipping Weight	25 g	Net Weight	24 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Length	18 cm		
Tipdimensions	1.4 x 1 mm		
Tips	Serrated		
Tipshape	Straight		