



Micro-Adson Forceps - 1x2 Teeth/Straight/12cm



General Info

Product Number	11019-12
Category	Forceps & Tweezers > Toothed
EAN	4251909501433



Find 11019-12 on
finescience.com

Technical Info

Shipping Weight	20 g	Net Weight	18 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Length	12 cm		
Tipdimensions	0.8 x 1.3 mm		
Tipthickness	1.3 mm		
Tipwidth	0.8 mm		
Tips	1x2 Teeth		
Tipshape	Straight		