



Surgical Standard Forceps - 1x2 Teeth/Straight/12cm



General Info

|                |                              |
|----------------|------------------------------|
| Product Number | 11021-12                     |
| Category       | Forceps & Tweezers > Toothed |
| EAN            | 4061706024109                |



Find 11021-12 on  
[finescience.com](https://www.finescience.com)

Technical Info

|                        |                                 |                 |      |
|------------------------|---------------------------------|-----------------|------|
| Shipping Weight        | 20 g                            | Net Weight      | 19 g |
| Package Size           | 1 Pc.                           | Unit of Measure | EA   |
| Minimum Order Quantity | 1                               |                 |      |
| Alloymaterial          | Surgical Stainless Steel (Inox) |                 |      |
| Length                 | 12 cm                           |                 |      |
| Tipdimensions          | 1.7 x 1.2 mm                    |                 |      |
| Tips                   | 1x2 Teeth / Standard            |                 |      |
| Tipshape               | Straight                        |                 |      |