



Surgical Standard Forceps - 1x2 Teeth/Straight/15.5cm



General Info

Product Number	11021-15
Category	Forceps & Tweezers > Toothed
EAN	4061706024147



Find 11021-15 on
finescience.com

Technical Info

Shipping Weight	35 g	Net Weight	33 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Length	15.5 cm		
Tipdimensions	2.5 x 1.7 mm		
Tips	1x2 Teeth / Standard		
Tipshape	Straight		