

F · S · T

## Surgical Narrow Forceps - Slim1x2 Teeth/Straight/14.5cm



## Allgemein

Produktnummer	11023-14
Kategorie	Forceps & Tweezers > Toothed
EAN	4251909501495



Finde 11023-14 auf  
[finescience.com](https://finescience.com)

## Technische Informationen

Versandgewicht	25 g	Nettogewicht	21 g
Paketgröße	1 Pc.	Maßeinheit	EA
Mindestbestellmenge	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Length	14.5 cm		
Tipdimensions	2 x 2.2 mm		
Tips	1x2 Teeth / Narrow		
Tipshape	Straight		